

WAK-Chemie Medical GmbH introduces its Ease-of-Rinse Cleaning Solution to OEM Community:

SystemClean SpeedFlush

Non-Foaming Acidic Cleaning Agent for
Fixed-Tip Liquid Transfer Robotics and Analyzers on a Daily Basis

SystemClean SpeedFlush achieves to provide excellent cleaning results by its proprietary formulation and easily removes of all residues from the system.

An organic acid solution combined with organic polymers of optimised chain length and surfactants overcomes the rinsing problems of analytical and transfer systems and excels in:

Cleaning Power in Cleaning Cycle or Overnight Incubation:

- detachment of sludge and other deposits from the interior of the system
- complete wash-away of cellular and other particles

Ease-of-Rinse:

Complete removal of all traces of the cleaning solution is achieved in 1-2 rinse cycles only, depending on the features of the liquid transfer system. Benefits include:

- elimination of carry-over effect
- drastically reduced downtime for daily maintenance routines

Highest Quality Ingredients:

- All ingredients are of pure grade and defined composition

Services for OEM Customers:

- WAK-Chemie Medical GmbH provides free samples for evaluation purposes
- comprehensive technical and scientific support for adaption and regulatory affairs



REF WAK-SCSF-500

500 ml Concentrate for 5 Liters Working Solution

L-01-158-02-0-00

Siemensstraße 9
D-61449 Steinbach/Ts. • Germany
Tel. +49 61 71 / 28 43 0 • Fax +49 61 71 / 69 56 46
e-mail: info@wak-chemie.com
www.wak-chemie.com



**WAK - Chemie
Medical GmbH**

